

SlideRightQA

AI-powered quality assurance for digital pathology



- **Automated detection of slide scan quality including edge cases**
- **Centralized dashboard with intuitive interface**
- **Vendor-neutral software solution designed for pathology QC staff**
- **Streamlined workflow efficiency**
- **Reduced turnaround times**
- **Reduced manual review hours**
- **Resource optimization**
- **Minimized rework expenses**

Cut costs, not corners in your digital pathology workflow

SlideRightQA supports continuous process improvement in your digital pathology lab by automating quality review of whole slide images.

The AI-powered software identifies issues such as out of focus and missed tissue, ensuring that only the highest quality images make it to your pathologists' queue. This helps you improve turnaround times as well as reduce costs and delays.

By eliminating the need for manual quality and compliance review of every whole slide image, SlideRightQA helps you streamline your workflows, increase efficiency and save time and money!