

Coronis 2MP (MFGD-2320)

2 MegaPixel premium diagnostic display system



Product specifications

CORONIS 2MP (MFGD-2320)

General specifications

Analog	VGA D-SUB15 connector, BNC
Approvals	Safety: <ul style="list-style-type: none">CE, UL60950, CSA C22.2 No 60950 (c-UL), CB, IEC 60950, DEMKO EN 60950 Emission: <ul style="list-style-type: none">EN 55022, EN 50082-1, FDA 510 k
Color imaging	No
Dark room contrast	1000:1 typical
Display technology	<ul style="list-style-type: none">TFT AMLCD VA Mode LCDGrayscale, 3 subpixels per pixel
DVI	Complying to VESA DVI
Ergonomics	<ul style="list-style-type: none">Tilt & swivel base with P/L rotation100 mm VESA mounting standardAR coated front cover
Gray imaging	True
Luminance	400 cd/m ² (116.7 fL) calibrated, 600 cd/m ² (175.1 fL)
Resolution	2MP (1600 x 1200 / 1200 x 1600)
Stabilization, calibration and automated QA	Embedded I-Guard
Supported resolution	<ul style="list-style-type: none">1600 x 1200 landscape1200 x 1600 portrait (Optional portrait Accelerator)VGA boot
USB	USB-hub with 1 up-and 2 downstream ports
Viewing angle	170° (at 10:1 contrast)

Last updated: 09 Jul 2024

© 2024 Barco nv. All rights reserved. Reproduction in whole or in part without written permission is prohibited. All brand names and product names are trademarks, registered trademarks or tradenames of their respective holders. Due to continued innovation, information and technical specifications are subject to change without prior notice. Please check www.barco.com for the latest specifications.