


# MDRC-2124

24-inch clinical review display



## **Consistent, calibrated images at the point of care**

The MDRC-2124's high-quality LCD panel provides high brightness, rich contrast and accurate display of colors and grayscales. The wide 24-inch screen creates a flexible workspace and is ideally suited for simultaneous use of multiple applications and palettes.

Advanced Backlight Output Stabilization (BLOS) technology ensures stabilized luminance over time.

Thanks to the DICOM Part 14 Look-Up Table, the MDRC-2124 is DICOM-compliant "out-of-the-box".

## **User-friendly Quality Assurance**

The MDRC-2124 is shipped with Medical QAWeb software for accurate Quality Assurance and remote display management (hardware calibration can be performed by means of an external sensor).

DynamicPDF