

# Coronis Fusion 4MP DL (MDC C-4130)

4 MegaPixel wide-screen pathology display system

- **Ultimate reading flexibility**
- **I-Guard technology for accurate images**
- **Automated QA & calibration**

Coronis Fusion 4MP DL is a 30-inch pathology display system designed to further enhance flexibility and productivity in pathology imaging.

### **Enhanced viewing comfort**

With the Coronis Fusion 4MP DL, there are no more boundaries. Its ingenious design allows you to use the 30-inch workstation as two seamless 2 MegaPixel heads or one wide-screen 4 MegaPixel display. The choice is yours.

Addressing the needs of both grayscale and color imaging, the Coronis Fusion 4MP DL allows you to read multiple images side by side on a single screen.

### **A perfect view from every angle**

Coronis Fusion 4MP DL features the latest, breakthrough In-Plane Switching (IPS-Pro) LCD technology. This technology brings the most advanced LCD viewing characteristics on the market, setting new standards for brightness and contrast, even from a wide viewing angle.

DynamicPDF