

Custom medical solutions

Designed with your needs in mind



- All-inclusive OEM service
- Individual guidance and follow-up
- Products delivered on time, on spec, and on budget

Everything you need. And more.

We develop:

- X-ray displays
- CT displays
- MRI displays
- Nuclear medicine displays
- Electrophysiology displays
- Ultrasound displays
- C-arm displays
- Interventional imaging displays
- Surgical displays
- Surgical robotics displays
- Patient monitoring displays
- User interface devices
- All-in-one devices

Devices that just feel right

Our user interface devices present information on an easy-to-use touch screen or smart display. Every device is engineered to meet the strictest medical



requirements and follows the latest design specifications. Easy to use, extremely flexible, and highly durable, they offer the best ergonomic user experience.

Built according to the most stringent protocols, these devices allow for rapid prototyping and validation. Expect all the capability you wish for in a device that just feels right.

Features of Barco's custom-built solutions:

- Outstanding performance (with patented breakthrough innovations)
- Perfect form-fit-function
- Proven solutions (tested in-house)
- Fully medically certified
- Effective time to market (in-house production line)
- Sustainable design
- On time and on spec