

RetroData 2100



Product specifications**RETRODATA 2100****General specifications**

Light output	635 Nit/185 fti
Resolution	VGA, SVGA (800x600), XGA (1024x768)
Contrast	>300:1 (full field contrast)

Last updated: 09 Jul 2024

© 2024 Barco nv. All rights reserved. Reproduction in whole or in part without written permission is prohibited. All brand names and product names are trademarks, registered trademarks or tradenames of their respective holders. Due to continued innovation, information and technical specifications are subject to change without prior notice. Please check www.barco.com for the latest specifications.