

# CliniScape (MTCP-0110)

10-inch ruggedized mobile clinical assistant



APPROVED



TOUCH



COLOR



10"  
SCREEN



**Product specifications****CLINISCAPE (MTCP-0110)****General specifications**

CPU	Intel® Core™ Solo
CPU Speed	1.2 GHz
RAM Size	2 GB DDR2 SDRAM
Battery capacity	3900 mAh
Hard disk	80 GB
Operating system	Microsoft® Genuine Windows® operating system
Technology	Transmissive TFT LCD
Active screen diagonal	264 mm (10.4")
Active screen width	210.4 mm (8.28")
Active screen height	157.8 mm (6.21")
Color imaging	Yes
Gray imaging	No
Resolution	XGA (1024 x 768)
Luminance	160 cd/m2 (typical)
Touch screen	Pen-sensitive digitizer type
Contrast	250:1 (typical)
Dimensions	265 x 250 x 29.5 mm (10.4 x 9.8 x 1.2") (W x H x D) Depth including handle: 78 mm (3.1")
Weight	1.8 Kg (4 lbs) Approx.
Approvals	EN / IEC 60950-1, UL 60950-1 EN / IEC 60601-1 CAN/CSA-C22.2 No.691.1-M90 FCC CFR 47 P.15 Sub B -Class B EN 55022 / CISPR 22 -Class B EN 55024 EN / IEC 60601-1-2 EN 55011/ CISPR11 -Class B FCC CFR 47 P.15 Sub C RSS-210, SAR, IC ID C-UL-US, CE, DEMKO, FCC
Operating temperature	0 to +35°C (32 to 95°F)
Ambient Light Compensation (ALC)	No
Uniformity Technology	No
Optional configurations	2D Barcode imager
Options	<ul style="list-style-type: none"><li>• Docking station</li><li>• VESA mounting adapter</li><li>• Stylus pen kit</li><li>• Spare battery</li><li>• Spare fan cartridge</li></ul>

**Last updated: 08 Jul 2024**

© 2024 Barco nv. All rights reserved. Reproduction in whole or in part without written permission is prohibited. All brand names and product names are trademarks, registered trademarks or tradenames of their respective holders. Due to continued innovation, information and technical specifications are subject to change without prior notice. Please check [www.barco.com](http://www.barco.com) for the latest specifications.