

Coronis Fusion 4MP DL (MDCC-4130)

4 MegaPixel wide-screen pathology display system



- **Ultimate reading flexibility**
- **I-Guard technology for accurate images**
- **Automated QA & calibration**

Coronis Fusion 4MP DL is a 30-inch pathology display system designed to further enhance flexibility and productivity in pathology imaging.

Enhanced viewing comfort

With the Coronis Fusion 4MP DL, there are no more boundaries. Its ingenious design allows you to use the 30-inch workstation as two seamless 2 MegaPixel heads or one wide-screen 4 MegaPixel display. The choice is yours.

Addressing the needs of both grayscale and color imaging, the Coronis Fusion 4MP DL allows you to read multiple images side by side on a single screen.

A perfect view from every angle

Coronis Fusion 4MP DL features the latest, breakthrough In-Plane Switching (IPS-Pro) LCD technology. This technology brings the most advanced LCD viewing characteristics on the market, setting new standards for brightness and contrast, even from a wide viewing angle.