

AC-WB02

Multimedia wallbox for Smart Terminals



Ideal for patient entertainment

The multimedia wall box increases the connectivity options of the Smart Terminal by providing an Ethernet connection at the bedside allowing patients and hospital staff to use their own PC notebooks and tablets.

The wallbox also potentially allows games consoles to be connected and even for patients to access their own music and photo files.

From a clinical perspective, hand held barcode scanners can be connected to the wall box for medicine dispensing and patient identification.

The unique ceiling mount option provides installers with optional mounting for the Smart Terminals in areas where space is limited or unsuitable for wallbox fixing.

A dual arm ceiling mount is also available and provides multiple installations of arms in confined areas.

Product specifications**AC-WB02****General specifications**

Height	Ceiling mounted fixture
Cover	Mounting plate covers
Rotation	Rotational stops to prevent cable damage
Cable management	Flexible cable installation and management Integrated cable covers for easy cable installation Covered cable joints and articulation points Large inner cross-section for cables
Color	Cool Gray 3C
Reach	Custom mounting lengths (on request)
VESA	Integrates with Arcos arms and flat screen monitors with VESA-standard
Movement	Perfect parallel guidance during horizontal movement Easy adjustment without dismantling the arm Easy arm movement in all articulations

Last updated: 08 Jul 2024

© 2024 Barco nv. All rights reserved. Reproduction in whole or in part without written permission is prohibited. All brand names and product names are trademarks, registered trademarks or tradenames of their respective holders. Due to continued innovation, information and technical specifications are subject to change without prior notice. Please check www.barco.com for the latest specifications.