


QAWeb Private Practice

Automated Quality Assurance of medical displays in private practices



QAWeb Private Practice guarantees optimum and stabilized image quality in every private practice. The software comes out-of-the-box (free of charge) with Barco's Eonis displays.

Ease of use

While IT administrators manage Quality Assurance for large healthcare enterprises, QAWeb Private Practice has been specifically designed for healthcare environments which do not have these dedicated resources. The software is highly intuitive and allows for automatic software updates in order to address your specific needs.

Fit for purpose

To ensure the best image quality, QAWeb Private Practice automatically suggests personalized settings for your display based on your user profile (i.e. occupation and geographical location). The software also features an easy-to-use tool to optimize image quality in changing ambient light conditions.

Consistent image quality

The front-of-screen consistency sensor on Barco's Eonis displays works seamlessly with QAWeb Private Practice to ensure a consistent image over time. It automatically stabilizes the image from the moment you switch on the display. What's more, QAWeb Private Practice provides you with instant feedback on the status of the display.

Why you'll love Barco's QAWeb Private Practice:

- Specifically designed for healthcare professionals
- Automatic updates
- Wizard to easily personalize your display
- Automated and online Quality Assurance
- Instant feedback on display status