

Coronis 2MP (MFGD-2320)

2 MegaPixel premium diagnostic display system



技术规格

CORONIS 2MP (MFGD-2320)

一般规格

模拟	VGA D-SUB15 connector, BNC
认证	Safety: <ul style="list-style-type: none">CE, UL60950, CSA C22.2 No 60950 (c-UL), CB, IEC 60950, DEMKO EN 60950 Emission: <ul style="list-style-type: none">EN 55022, EN 50082-1, FDA 510 k
彩色成像	No
暗室对比度	1000:1 typical
显示技术	<ul style="list-style-type: none">TFT AMLCD VA Mode LCDGrayscale, 3 subpixels per pixel
DVI	Complying to VESA DVI
人体工程学	<ul style="list-style-type: none">Tilt & swivel base with P/L rotation100 mm VESA mounting standardAR coated front cover
灰色成像	True
亮度	400 cd/m ² (116.7 fL) calibrated, 600 cd/m ² (175.1 fL)
分辨率	2MP (1600 x 1200 / 1200 x 1600)
稳定性, 校正和自动化质量控制	Embedded I-Guard
支持分辨率	<ul style="list-style-type: none">1600 x 1200 landscape1200 x 1600 portrait (Optional portrait Accelerator)VGA boot
USB	USB-hub with 1 up-and 2 downstream ports
视角	170° (at 10:1 contrast)

生成于: 09 Jul 2024

© 2024 巴可公司。保留所有权利。未经书面许可, 禁止全部或部分复制。所有品牌名称和产品名称均为其各自持有者的商标、注册商标或商号。由于不断创新, 信息和技术规格如有变更, 恕不另行通知。请查看www.barco.com.cn以了解最新规格。