


# MD-4221

42-inch High Definition LCD display for the operating room



## Bringing PACS to the OR

Fitted with a 10-bit high-contrast LCD panel and a built-in DICOM Look-Up Table, the MD-4221 is optimal for consulting PACS images in the OR.

## User-definable screen layout

The MD-4221 lets you combine two images - in any mix of data and video - next to each other or as Picture-In-Picture. Its wide viewing angle of 178° ensures that all this information is clearly visible for everyone in the surgery suite.

## Designed for the operating suite

An anti-reflective protective front glass allows the monitor to be cleaned and disinfected (IPX1 compliant), while its fan-less design avoids the spread of contaminants in the operating room.