

# Coronis Fusion 4MP (MDCC-4330)

4MP pathology display system

- **Color fusion display for flexible pathology imaging**
- **With toolset to boost productivity**
- **Automated QA and calibration**

## **A pleasure for your eyes**

Coronis Fusion 4MP features an exclusive LCD panel with the latest IPS technology, offering best-in-class image quality. Its high brightness and high contrast ratio help unveil the subtlest of details. Its large 30" screen real estate, wide aspect ratio, and full ACR compliance make this display perfect for pathology imaging.

## **Boost your productivity**

Coronis Fusion 4MP comes with smart features that enhance control and productivity. Unique SpotView™ technology allows you to boost the luminance in a region of interest while dimming the surrounding image area so you can focus on more details more efficiently.

Like all of Barco's medical display systems, Coronis Fusion 4MP comes with MediCal QAWeb, a cloud-based technology for automated calibration and Quality Assurance to ensure maximum uptime of the display with no need for human intervention.

## **Designed with ergonomics in mind**

With Barco's Coronis Fusion 4MP, switching between workstations becomes a thing of the past. You can use the single-head workstation as two bezel-free heads or as one wide-screen workspace. The choice is yours.

Its widescreen, ergonomic form factor, wide viewing angle, and DimView™ feature (to automatically dim auxiliary displays) have been designed for the way you work and help you focus on the smallest image details.

## **A one-stop-shop solution**

Barco medical displays include proprietary display controllers that are validated with the latest workstations and with all major PACS applications. All components are under the full 5-year warranty (including the display backlight) for complete peace of mind.