


# MFCD-1219 TS

19-inch medical color display with touch screen interface



© 2008 BARCO Medical Systems Group

**Product specifications****MFCD-1219 TS****General specifications**

Technology	AMLCD dual domain IPS
Active screen diagonal	483 mm (19")
Active screen width	376.3 mm (14.8")
Active screen height	301.1 mm (11.9")
Color imaging	Yes
Gray imaging	No
Resolution	1.3MP (1280 x 1024)
Luminance	220 cd/m <sup>2</sup> (Typ.)
Contrast	800:1 (Typ., dark room)
Viewing angle	178° (Hor./Vert.)
Video inputs	D-SUB15 HD, DVI
Dimensions	450 x 463 x 214 mm (17.7 x 18.2 x 8.4") (W x H x D, unpacked)
Weight	8.4 kg (18.5 lbs)
Approvals	CE, UL 60950, c-UL, IEC 60950, FCC-B, iCES 003, EN 55022 limit B, CE, EN 55024
Operating temperature	0 to +40°C (32 to 104°F)
Ambient Light Compensation (ALC)	No
Uniformity Technology	No

Last updated: 08 Jul 2024

© 2024 Barco nv. All rights reserved. Reproduction in whole or in part without written permission is prohibited. All brand names and product names are trademarks, registered trademarks or tradenames of their respective holders. Due to continued innovation, information and technical specifications are subject to change without prior notice. Please check [www.barco.com](http://www.barco.com) for the latest specifications.