

TV & radio (JAO SW-CL-TV^R)

Bedside computing software for patient entertainment

A better patient experience

With Barco's bedside computing software, patients can enjoy their favorite TV or radio show while they're in the hospital. What a great way to enhance their hospital experience! With the TV & radio module, any TV or radio channel can be distributed via the IP TV head end.

These programs can be received freely or encoded. The Barco solution includes the option of central decoding of encoded programs. TV stations can be provided free of charge or for a fee. Packages are also possible, such as some channels free and others chargeable.

Barco's TV & radio software offers:

- TV (preview window, full screen mode, volume control...)
- Radio (volume control, channel previews...)
- TV/radio billing (overview of purchased packages)
- Hardware buttons (channel up/down, volume control...)

Barco's bedside computing software is ideal for use on its JAO Smart Terminal range and runs on a wide variety of hardware and third party systems.