

# MR-100-E

MR-100-E mission recorder



**Product specifications****MR-100-E****General specifications**

Storage capacity	4x500 GB removable HDD max
Video Input	Dual Analog & Digital Video Input • 800x600 to 1920x1200 • MPEG-4 encoding
Video Output	Safe pass-through • Dual Analog & Digital Video output • 800x600 to 1920x1200
Audio Input	Analog Stereo Audio Input, 600 $\square$ Sampling 16-bits 44.1 KHz
Audio output	Analog Stereo Audio Output (2W/ 8 $\square$ )
Interfaces	USB 2.0, Gb Ethernet
Operating system	Windows XP/XPe/7, Linux
Power supply	115/230 VAC auto-range
Dimensions	483 x 89 x 482mm (with front handles)
Weight	< 14 kg
Vibration	MIL-STD-167 and MIL-STD 810F
Humidity	MIL-STD-167 and MIL-STD 810F
Operating temperature	0°C, +50°C according to MIL-STD 810F
EMI / EMC	EMI/EMC CE/FCC class A

Last updated: 09 Jul 2024

© 2024 Barco nv. All rights reserved. Reproduction in whole or in part without written permission is prohibited. All brand names and product names are trademarks, registered trademarks or tradenames of their respective holders. Due to continued innovation, information and technical specifications are subject to change without prior notice. Please check [www.barco.com](http://www.barco.com) for the latest specifications.